



## Serial purification of phages by streaking

(Adapted from Kauffman and Polz 2018)

### Introduction

After isolating a phage from a plate and before producing a high titre stock, it is necessary to re-isolate it at least twice.



### Materials (example for *V. crassostreae*)

- Marine broth (MB)
- Marine agar (MA)
- Top agar (entre 0.4 et 0.2% selon les phages)

### Protocol:

1. Boil heavily the top agar bottle and check that the media is completely liquid and homogeneous. Keep at 55°C with stir.
2. Dispense 100 µl of host overnight culture on dry MA plates with a 2,5 ml fresh agar overlay (if top agar has "clumps" discard it or boil again). **Next steps have to be done quickly since top agar has to be molten.**
3. Quickly, insert yellow tip into the plaque or liquid phage source and spread it (from the border of the plate to the centre) on the still-molten agar overlay (blue line on figure).
4. Use a second yellow tip to make three separate strokes (Z) from the initial zone through the top agar (**never come back when you are doing the Z!**) (red line on figure).
5. Finally, use a third yellow tip to make a new Z through the strokes done in the step 4 (green line on figure).
6. Incubate the plates for 24h. Once plaques appear, use this plate as phage source and repeat the process from step 2 twice.
7. Collect a single plaque using 1 ml tip (previously, cut a little bit the 1 ml tip) from the third serial purification. Eject the agar plugs of plaque in Eppendorf with 750 µl of MB and allow the elution of phage particles into the media overnight at 4°C.

8. After soaking, use a syringe to filter the phages by 0,2  $\mu\text{m}$  in a new tube (you will lose volume due to the filtration, approx. your final volume will be 300  $\mu\text{l}$ ). Purified phage stock could be stored at 4°C (or -80°C with 25% glycerol) or use for preparation of high titer stocks.

**Important tip!**

Most of the phages in our collection (Caudovirale) are isolated using 0.4% TOP agar. Schizotequatroviruses require 0.2% TOP agar made with VWR agar J637-500G.

**Reference**

Kauffman KM, Polz MF. Streamlining standard bacteriophage methods for higher throughput. *MethodsX*. 2018 Jan 31;5:159-172. doi: 10.1016/j.mex.2018.01.007.